Twin Art

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COMPOSITION

Polyethylene frame padding LDPE, nylon buttons, closure with velcro straps.

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INTENDED USE

Immobilization of the cervical spine, minimizes flex extension.

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WARNING

Don't use collars of inadequate measures.

Always refer to your doctor before use and to the instructions on the use of medical device

100% Radiolucent - non sterile - single patient use - Latex Free - Single use

"It is advisable to report to Artsanity Srl and the competent national authority any serious accident occurring in relation to the device"

Product manufactured and distributed in accordance with Regulatory (UE) 2017 / 745

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Dispositivo monouso

STORAGE AND CLEANING

Store in a cool, dry, clean place away from direct sunlight and heat sources.

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Separate two elements. Take the front and insert the plugs.

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Stabilize the patient's head. Pull off the collar to the chin and to the sternum, keep it in place and wrap the neck.

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Girding neck with band and secure it with velcro while keeping neck position.

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Place the back behind head centered manner and bring together in velcro band.

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
Tighten velcro band to stabilize the patient's neck position.

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Artsanity srl

Via della Borsa, 18 – 31033 Castelfranco Veneto (TV) – Italy - Tel. 0039.0423.494209 – Fax 0039.0423.495866

C.F. e P.I. 03190310262 - Cap. Soc. € 10.330,00 - Reg. Imp. di Treviso nr. 03190310262 - R.E.A. n. 227049

www.artsanity.it – e-mail: info@artsanity.it – PEC: artsanity@legalmail.it

Twin Art

Polietilene bianco "citotoxicity Test According ISO 10993-5"

TEST DI CITOTOSSICITA' DEL MATERIALE Rif. NORMATIVA ISO 10993-5
MATERIAL CITOTOXICITY TEST ACCORDING ISO 10993-5

Test report number: ISO 202002-02852I_engl

Test Report

*Toxicity testing of a material extract,
Cytotoxicity test*

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Test material received: Mar 09th 20

Test performed: Apr 17th 20

Result: The polietilene bianco altezza mm. 6, Code LD6BI, LOT9
did not cause a cytotoxic effect.

CYTOX biological testing of medical devices

Test report: toxicity testing of a material extract Apr 29th 20

page 1/3

Artsanifty srl

Via della Borsa, 18 – 31033 Castelfranco Veneto (TV) – Italy - Tel. 0039.0423.494209 – Fax 0039.0423.495866
C.F. e P.I. 03190310262 - Cap. Soc. € 10.330,00 - Reg. Imp. di Treviso nr. 03190310262 - R.E.A. n. 227049
www.artsanifty.it – e-mail: info@artsanifty.it – PEC: artsanifty@legalmail.it

Test report number: ISO 202002-02852I_engl


Description of the test procedure:

Normative References: ISO 10993-1 (2018); ISO 10993-5 (2009); ISO 10993-12 (2012)

The material was extracted for 24 h at 37°C and a partial pressure of 5 % carbon dioxide in extraction medium (DMEM medium with antibiotics, without fetal calf serum [FCS]). The material / extraction volume ratio was 0,2 g material per ml extraction medium. After extraction the extraction medium was sterile filtered and supplemented with sterile FCS (concentration of FCS in extraction medium: 10 %). The FCS-supplemented extraction medium was pipetted under sterile conditions on precultivated cells of the mouse fibroblastic cell line L929 and incubated for 48 h at 37°C and a partial pressure of 5 % carbon dioxide. The extract was tested in four dilutions (90 %, 30 %, 10 % and 3,3 %). Each dilution was tested in four parallel experiments.

Triton X 100 was used as a toxic control substance (concentration in the experiment: 1 % v/v). Cell culture medium was used as a non-toxic control. After the 48 h incubation period the protein content of the cell culture was determined by the method according to Bradford.

Results:

Result data	Protein content n=4					
µg/ml	Triton 1 %	Culture medium	Extract 90 %	Extract 30 %	Extract 10 %	Extract 3,3 %
Mean	2,16	48,90	48,36	48,55	48,51	49,44
Standarddev.	1,23	1,66	0,30	1,51	1,25	0,38

In the presence of Triton X 100 in the cell culture medium 4,4 % of the protein content compared to the negative control was reached. This value is within the valid range of 15 % protein content or less compared to the negative control.

CYTOX biological testing of medical devices

Test report: toxicity testing of a material extract Apr 29th 20

page 2/3

Artsanity srl

Via della Borsa, 18 – 31033 Castelfranco Veneto (TV) – Italy - Tel. 0039.0423.494209 – Fax 0039.0423.495866
C.F. e P.I. 03190310262 - Cap. Soc. € 10.330,00 - Reg. Imp. di Treviso nr. 03190310262 - R.E.A. n. 227049
www.artsanity.it – e-mail: info@artsanity.it – PEC: artsanity@legalmail.it

Twin Art

Polietilene bianco "citotoxicity Test According ISO 10993-5"

Test report number: ISO 202002-02852I_engl

Materials are considered cytotoxic, if the material extract leads to a protein content of the test cells of less than 70 % compared to the negative control. This was not the case in this test. The material extract did not show a cytotoxic effect.

Result: The polietilene bianco altezza mm. 6, Code LD6BI, LOT9
did not cause a cytotoxic effect.

Explanatory notes:

none

Test performed by: Dr. D. Scheddin

authorized by: Dr. D. Scheddin
(Dr. D. Scheddin / CEO CYTOX)

It is not allowed to publish only parts of this test report without written approval of CYTOX.

CYTOX biological testing of medical devices

Test report: toxicity testing of a material extract Apr 29th 20

page 3/3

Artsanifty srl

Via della Borsa, 18 – 31033 Castelfranco Veneto (TV) – Italy - Tel. 0039.0423.494209 – Fax 0039.0423.495866
C.F. e P.I. 03190310262 - Cap. Soc. € 10.330,00 - Reg. Imp. di Treviso nr. 03190310262 - R.E.A. n. 227049
www.artsanifty.it – e-mail: info@artsanifty.it – PEC: artsanifty@legalmail.it